
H-Max S
H-Max M
H-Max C

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**PRODUCT DESCRIPTION
& SURGICAL TECHNIQUE**

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The H-MAX stems develop the principles of the Corail Philosophy; the stems are straight with a rectangular cross sectional area and HA coating for a fast and effective osteointegration. The system consists of three stem types (Fig. 1):

H-MAX S

Monoblock stem in the version standard (CCD 134°) and lateralized (CCD 131°).

H-MAX M

Modular stem which can be coupled with one of the 12 modular necks, allowing to choose the implant features in terms of:

- Offset
- Length
- Anteversion

H-MAX C

Cemented monoblock stem in the version standard (CCD 134°) and lateralized (CCD 131°).

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Materials

A specific study to ensure primary and secondary stability led to the following choices in terms of materials and coatings.

H-MAX S / H-MAX M

- The stems are made of forged and heat treated Ti6Al4V titanium alloy.
- The external surface has a macro-roughened surface that is obtained by blasting with a controlled flow of corundum particles. By controlling pressure, speed and grain size of the jet it is possible to create an uniform surface roughness ($Ra\ 3\div 6\ \mu m$) in order to enhance the bone integration. To ensure a fast bone growth around the prosthesis, a layer of hydroxyapatite ($55\ \mu m$) is applied along the whole length of the stem. Necks for H-MAX M are made in CoCrMo.

H-MAX C

The stem is made in M30 (FeCrNiMnMoNbN) and mirror-polished.

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[4]

► Fit and fixation in the meta-epiphyseal area, thanks to a wide medial curvature characterised by a 100 mm radius that guarantees an excellent anatomical fit, which is indispensable for primary and secondary stability (Fig. 2).

► A-P profile, characterised by the "V" shape in the proximal part with an 8° angle that guarantees primary stability of the implant (Fig. 3). This design provides a proximal distribution of the torsional strains, supporting flexional stresses through diaphysis.

► Maximum torsional stability achieved thanks to the rectangular section and a marked lateral flare.

► Tapered profile of the distal part to improve adaptability to the most common bone morphologies and to facilitate insertion of the stem in the canal, avoiding irregular transfer of the load to the bone structure. This minimizes the undesired stress shielding and the point effect that is responsible for thigh pain (Fig. 4).

The design of the stem, self locking, is optimized in the proximal portion to increase the mechanical stability and ensure an adequate distribution of the loads on the bone surface.

For the uncemented stems H-MAX M and H-MAX S: the horizontal ridges in the proximal part of the stem avoid sinkage and enhance the secondary stability of the implant.

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[5]

The distal region is characterized by vertical ridges which increase the rotational stability of the stem (Fig. 5).

► All necks are lowered and tapered to avoid accidental contact between stem and acetabular cup, increasing the range of motion. Modular ceramic or metal heads of various lengths and diameters can be used on the 12/14 taper.

► The necks are mirror-polished to reduce abrasion of the polyethylene cup in the event of abnormal contact.

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
[6]

Different combinations
of modular necks

- a. Standard - Lateralizing
b. Anteverted - Standard - Retroverted.
c. Short - Long.

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[7]

In the H-MAX M stem the modularity has been achieved with a system of modular necks (Fig. 6) which guarantees:

- Less variations in cup positioning or stem anteversion;
- Reduced risk of dislocation and impingement;
- Improved adaptability of the implant to the different anatomic-physiological characteristics of the patient.

The 12 necks are the result of the possible combinations in terms of Length, Offset and Anteversion (Fig. 7).

Length: The necks are available in the versions Short and Long; which differ of 10.5 mm length along the axis of the neck. Using head sizes S,M, L it is possible to increase the length linearly, with 3.5 mm step, restoring the right tissue tensioning and the centre of rotation.

Offset: The necks are available in the version Standard and Lateralizing. The Lateralizing version is obtained increasing the offset of 5 mm.

Anteversion: The necks are available in Neutral, Anteverted (+10°) and Retroverted (-10°) versions.

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Indications and contraindications

Indications

Cementless and cemented straight stems for partial or total hip replacement in case of:

- Osteoarthritis
- Traumatic arthritis
- Rheumatoid arthritis
- Avascular necrosis of the femoral head
- Traumatic fracture of the femoral neck

Contraindications

- Acute or chronic infections, local or systemic infections
- Serious muscular, neurological or vascular disease affecting the concerned limb
- Bone destruction or poor bone quality, which could compromise implant stability
- Any concomitant disease that might affect the implanted prosthesis
- Allergy to material
- In case of Metal-on-metal couplings: patients with renal impairment
- BMI > 28 (for H-MAX M Stem)

Warnings

With the “long version” modular necks the stem H-Max M, only S, M & L lengths of femoral heads can be used (both ceramic and metal heads). It is not allowed to use XXL and XXXL heads with any of the modular necks.

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Conditions which can adversely affect the outcome of the operation:

- Severe osteoporosis
- Severe bone deformities
- Local bone tumours
- Systemic diseases and metabolic disorders
- Patient's history of infections or falls
- Drug addiction and/or drug, alcohol or medicine abuse
- Obesity
- Heavy physical activity associated with frequent shocks in which the prosthesis can be excessively overloaded

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
Sizes

The H-MAX M is manufactured in 10 different sizes. The CCD is 134° with the Standard necks.

H-MAX S and H-MAX C are manufactured in 11 sizes for the two versions Standard and Lateralizing. The CCD in the Standard version is 134°; the CCD in the Lateralizing version is 131°.

The dimension of the stems grows harmonically with size; in the frontal view the increase in width is 1 mm for each size (Fig. 8). In the lateral view the thickness increases of 0.5 mm for each size (Fig. 9).

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
[2]

1. Pre-operative Planning

To obtain the best results, pre-operative planning is recommended with the use of templates (showing a 15% enlarged image of the profiles). Use good quality frontal and axial view radiographs with adequate contrast that are large enough to contain the entire length of the pre-op template stems (Figs. 1-2).

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[4]

2. Choice of Stem Size

To find the right stem size the contour lines of the proximal part should fill the femoral epiphyseal region.

Make sure that the stem does not interfere with the femoral curvature along the axial projection.

It is important to consider that the fixation of the stem is obtained in the proximal part, therefore, the filling of the canal is not required in the distal part.


Identify the size and type of neck to restore the correct centre of rotation, checking the anteversion level in the sagittal plane.

Identify the neck resection level by making the apex of the greater trochanter coincide with the centre of the medium femoral head.

In case of implant of H-MAX S and C the templates indicate the centre of rotation in the two versions for the different lengths of the heads (Figs. 3-4).

For H-MAX M the templates indicate the centre of rotation for the different versions of the necks (Standard-Short, Standard-Long, Lateralizing-Short, Lateralizing-Long) and for the different heads lengths (Fig. 5).

For H-MAX M the templates also indicate the centre of rotation for the versions Standard-Short and Standard-Long in the axial projection (Fig. 6).

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
Note. Pre-operative planning provides useful indications for the use of the the implant but does not definitely extabish the size of the stem to be used, which must be verified during the surgery.

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[7]

3. Position of the Patient


Note. The first picture of the surgical technique shows the posterolateral approach with the patient in the lateral decubitus position. The surgical steps that follow also refer to the supine position and any other surgical access.

Patient in lateral decubitus position: select a posterolateral incision line (Fig. 7). Once the fascia lata is opened and the external rotators are detached, the joint capsule is cut and the head of the femur dislocated posteriorly. The top of the head is then exposed.

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[8]

4. Neck Resection

After having dislocated the femur, the femoral head is resected (Fig. 8).

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[9]

5. Reaming and Broaching

Start preparing the canal with the box chisel, which allows to open the great trochanter area (Fig. 9).

The canal is then opened using a reamer (Fig. 10).

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
The H-MAX instrument set consists of a series of broaches (Fig. 11); the broaches are connected to the broach handles. Broach handles are available in the following versions: "straight", "single offset" and "double offset" fitting to different surgical approaches (Fig. 12 a, b, c).

Connect the broach to the handle by lifting the lever and inserting the broach with its medial part oriented towards the lever (Fig. 13).

Tighten the lever until it is fully closed (Fig. 14).

Start broaching the canal using the smallest broach, maintaining the correct anteversion (approx. 15°).

Insert the broach until the tilted plane of the broach coincides with the resection line of the neck.

Note. Trial reduction with broach #08 can be done for H-MAX S and H-MAX C stems only (H-MAX M stem smallest size is #09).

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
Note. The resection line is defined during the Pre operative planning with the use of templates. If the resection line results to be different from the planning this difference has to be taken into account (Fig. 15 a, b, c).

Proceed with the next size up broach and continue until the epiphyseal seat is filled properly.


Once the best broach size has been achieved (which may not necessarily be the one planned during the pre-operative planning), remove the handle and leave the broach in situ (Fig. 15 b).

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6. Trial Reduction

Normally the acetabular time precedes the femoral time. Therefore at this point a trial reduction can be performed.

For H-MAX M stem the instrument set includes the trials of the 12 modular necks Short (red) and Long (grey).

For the H-MAX S e C the instrument set includes the trials for the necks in the two versions; Standard (green) and Lateralizing (blue) for each size. Insert the required trial neck with the neck positioner (Figs. 16-17).

Insert the trial head (Fig. 18) and perform the trial reduction (Fig. 19).

Remove the trial neck using the neck extractor.

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7. Definitive Stem Insertion

H-MAX M

Take the H-MAX M stem out of the box with the size corresponding to the last broach used and remove it from its sterile packaging (Fig. 20).

Using the broach handle, impact the stem, until the tilted plane of the stem reaches the level of the broach (Fig. 21c).

The neck can be further checked by inserting the trial neck on the final stem.

Clean and dry the conical taper of the stem.


Take the definitive neck from the sterile packaging and impact it with the neck positioner (Fig. 22) beating it on its axis.

[21a]

[21b]

[21c]

H-MAX C

After inserting the distal plug fill the femoral canal with acrylic cement from distal to proximal. Proceed with the introduction of the stem as for the H-MAX S as follows.

Please look at the LIMA-CMT 1, 2 and G brochure for indication of the distal plug positioning and the cementification of the canal.

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
H-MAX S

Take the H-MAX S stem out of the box with the size corresponding to the last broach used and remove it from its sterile packaging.

Screw the stem positioner on the definitive stem (Fig. 23).


The washer indicates the locking/unlocking direction (Fig. 24). The T Wrench, inserted in the stem positioner, is useful to lock up the stem to the stem positioner (Fig. 25). After removing the T wrench it is possible to insert the stem by beating on the stem positioner. The line of separation between the polished finishing and the HA coating corresponds to the level reached by the last broach used (Fig. 26).

Note. In the H-MAX C stem the wedge at the neck base corresponds to the level reached by the last broach used.

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

[30]

8. Final Head Insertion

The head length can be checked once again by using the trial heads.

Take a head with the determined diameter and length from its sterile packaging (Fig. 27).

Carefully clean and dry the neck taper, this is essential when using ceramic heads, fix the head by pressing it in and rotating it along the neck axis (Fig. 28). Remove the protective cap and tap slightly along the cone axis using the femoral head beater (available upon request) (Fig. 29). After having cleaned the articular surfaces, reduce the articulation (Fig. 30).

Note. For H-MAX M, with all “Long” version necks, the use of XL, XXL, XXXL head lengths is not allowed.

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[32]

9. Component Removal

If necessary, the prosthetic components can be removed.

To remove the femoral head, simply tap the bottom of the head in the axial direction using a beater.

Note. If only the head needs to be replaced, never use a ceramic head on the same cone.

H-MAX M

Remove the modular neck with the neck extractor (Fig. 31).

Lock the broach handle on the broach and beat backward to extract the stem.

H-MAX S / H-MAX C

Screw the stem positioner on the stem. Screw the stem extractor on the stem positioner and use the inertial beater to extract the stem (Fig. 32).

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Lima Orthopaedics New Zealand Ltd
Lima Orthopaedics UK Limited