SL Self Locking

Femoral Stems
Steli femorali

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Product Description

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The SL hip prosthesis relies upon the biomechanical concepts of the self-locking mechanism.

The straight stem has been flattened antero-posteriorly and primary stability is achieved by introducing the stem into the femoral shaft.

The introduction of the stem and the extrusion of cement in excess are favoured by the presence of two grooves on the anterior and posterior surfaces of the stem itself.

The cementless version is Plasma Spray coated with highly porous Titanium increasing osteointegration.

12/14 taper is compatible with metallic or ceramic heads for total implants and with endo or bipolar heads for applications in traumatology.

ADVANTAGES

- Simple and reliable system
- Possibility to choose the implantation technique during the operation
- Easy-to-use, precise instruments, similar for all versions

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
1. Preoperative Planning

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[1]

To obtain best results, preoperative planning using the given templates and x-rays is highly recommended. The profiles are enlarged by 15%, on the templates.

Use good quality frontal view and axial view x-rays with adequate contrast. The size of the x-ray must be large enough to contain the entire length of the stems on the templates.

The template must be placed according to the centre of the femoral head to obtain a precise neck resection level (Fig. 1).

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
2. Position of the Patient

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[2]

Note. The first picture of the surgical technique refers to the postero-lateral approach with patient lying on lateral decubitus. All the other surgical steps can equally refer to the supine position and to whatever surgical approach.


The patient lies on lateral position. Postero-lateral incision (Fig. 2). After opening the fascia lata and detaching the lateral rotators, the joint cap is cut and the head of the femur dislocated posteriorly. The top of the head is then exposed.

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3. Neck Resection

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[3]

Remove the femoral neck by trans-cervical resection, according to what planned preoperatively, at a 45° angle with the axis of the femur (Fig. 3).

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4. Preparation of the Medullary Canal

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Canal Chisel
Fustella per canale diafisario
(not supplied / non fornita) [4]

The diaphysis is punched along the femoral axis by means of the **canal chisel** (no supplied with instrument set) (Figs. 4-5).

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[5]

The preparation of the medullary canal starts using the smallest **rasp (A)** and keeping the correct anteversion (ca. 15°).

Open the hand lever of the **modular rasp handle (F)** and insert the rasp keeping the medial side towards the lever. Clamp the lever to the end (Figs. 6-7).

The modular rasps are used as trial prostheses.

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[7]

[8]

Complete the preparation of the prosthesis seating (Fig. 8).

Open the lever and remove the rasp handle leaving the rasp in the medullary canal (Fig. 9).

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[9]

5. Trial Reduction

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Keeping the rasp in situ, the determination of the head size is now possible.

Fit the support for trial heads SL and the trial head (C-E) (Figs. 10-11) over the rasp trunnion and reduce the joint seating the head into the acetabular cup previously implanted (Fig. 12). Choose the length of the femoral head to be used (S, M, L) testing the joint tension.

At the end remove the trial head, the support for trial head and the rasp by means of the handle.

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[13]

A special positioner -impactor (B) (Fig. 13) is available to fit into the medullary canal the stem of the selected size which is as the same as the latest rasp used.

To lock the stem on the positioner-impactor insert the taper of the stem into the adaptor terminal end and firmly screw the handle of the tool (Fig. 14) The tip gets in contact onto the proximal groove (Fig. 15).

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7. Cement Injection

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FOR CEMENTED STEMS ONLY:
Acrylic cement is filled into the medullary canal, after positioning the eventual distal plug (available on request) (Fig. 16).

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8. Stem Introduction

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Insert the stem into the medullary canal and gently tap on the head of the positioner –impactor until prosthesis is fully seated (Fig. 17).

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Note. For the cemented stems, manually remove the cement excess before the complete polymerization of the acrylic cement is achieved.

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[19]

[20]

Fit the trial femoral heads (no supplied with instrument set - Fig. 19) and reduce the joint seating the head into the acetabular cup previously implanted.

Choose the length of the femoral head to be used (S, M, L) testing the joint tension (Fig. 20).

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10. Head Insertion

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[21]

Take the definitive head in according to trialed one out of its sterile package.

Before to insert the head, please clean and wipe with care the neck taper. This precaution is very important when ceramic heads are used.

Fix the head on the neck pushing it along the neck axis and rotating it (Fig. 21).

Remove the protective shell (Fig. 22).

Gently beat with a plastic beater (no supplied with instrument set) along the taper axis (Figs. 23-24).

Note (according to EN ISO 21534):
The cemented SL stem can't be coupled with CoCrMo modular head.

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11. Final Reduction

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[25]

At the end check that coupling surfaces are cleaned and wiped before definitively reduce the joint (Fig. 25).

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12. Removal of the Components

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If necessary, all the components can be removed. For remove only the femoral head is enough beat with a punch (non supplied with instrument set) the head base along the axis.

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
Warning! Do never substitute on the taper a second ceramic head. In this case also the stem must be substituted.

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This image shows a full page of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page, typical of notebook or legal stationery. There are no margins, text, or other markings on the page.

[illegible]

For further information about our products, please visit our web site www.lima.it

This documentation has been written by Lima-Lto Sales & Technical Department

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Lima Austria GmbH
Lima Netherlands B.V.